Technical Data Sheet Revised:2020-05-08 | 1/2



**Art.-No.:** 6.7405.0401 6.7405.0402 6.7405.0404

**7405** 

# **Profi Primer Spray**

## Characteristics

Fast drying one-component self etching primer with active protection against corrosion.

Particularly suitable as an adhesion promoter for minor repairs and surfaces which have been partially sanded to the bare metal.

## Weldable!

Adhesion on iron, steel, zinc and aluminium (with proper pretreatment).

High protection against corrosion. Professional quality!

## V.O.C. Value

EC limiting value for the product (cat. B/e): 840 g/l (2007) This product contains a maximum of 666 g/l VOC [ 5.56 lbs/gal]

## Tech Tip

High covering power, fast drying, excellent adhesion. Recoatable with all commercial solvent and water based 1K- and 2K-Topcoats.

Weldable according to DVS-Info. 0501-No.275/670/99/925

Resistant to silicone remover. Active protection against corrosion.

Do not recoat with polyester spray fillers!



Technical Data Sheet Revised:2020-05-08 | 2/2



**Art.-No.:** 6.7405.0401 6.7405.0402 6.7405.0404

**7405** 

## Application

Substrate:

Iron, steel, zinc and aluminium

Do not process on thermoplastic surfaces.

Processing conditions: From +10°C and up to 80% relative air humidity.

Preparation:

Clean with 4CR silicone remover.

Shake can vigorously until the mixing balls move freely!

Application:

Spray to test. Spraying distance approx. 20 - 30 cm. 2 - 3 coats.

Flash-off:

3 - 5 min. between each coating

After use:

Turn can upside down and spray until the valve is clean, this prevents the valve from clogging up.

Further processing:

Dry sanding: P 400 - P500

## **Technical Data**

Colour red-brown (7405.0401) grey (7405.0402) or Black

(7405.0404)

 $\begin{array}{lll} \text{Dry-film thickness} & 50 - 60 \ \mu\text{m} \\ \text{Dust-dry} & 5 \ \text{min.} \\ \text{Dry to handle} & 10 \ \text{min.} \\ \text{Recoatable} & 30 \ \text{min.} \\ \text{Sandable} & 30 \ \text{min.} \\ \end{array}$ 

# Storage and Transportation

In tightly closed original containers at least 2 years shelf life.